




DreamStar™ Info & Auto



Quick user guide

Welcome

Congratulations for the purchase of your new DreamStar™ Info or Auto. Please refer to this quick user guide but always read the user manual prior using your DreamStar™ Info or Auto device.



Overview

The DreamStar Info & Auto can be delivered in 3 different configurations:

- Standard, with no heated humidification system on board.
- Evolve with cover, the humidification system is present but not active. Optional the heated humidifier can be activated and replace the cover by the water reservoir.
- Evolve with reservoir, the humidification system is active.

The DreamStar is delivered with the following parts:

DreamStar Info or Auto device, power supply cable, 22mm patient tubing, transport bag, data card, spare inlet filter and patient manual.

Optional

Fine filter, 15mm patient tubing, water reservoir on DreamStar Evolve with cover, cover on DreamStar Evolve with water reservoir.

Preparation for use



Setting up your DreamStar device

1. Place the device on the side of your bed, in a stable surface (eg bedside table).
2. Make sure that it is carefully and safely positioned.
3. Attach the cover or water reservoir onto the device.
4. Connect the power supply cable on the side of the DreamStar.
5. Plug the cable into the mains.
6. Connect the tubing onto the air outlet.
7. Connect the mask on the other part of the tubing.

Note:

- *Fill the water reservoir with water prior attaching it on to the device*
- *Respect the maximum water level indicator.*

WARNING

Make sure that the water reservoir is empty while traveling with your DreamStar Info or Auto device. There is a risk water entering the device if the water reservoir is not empty.

Setting up your DreamStar

Heated humidifier indication

Indicates that heated humidification feature is active and also show the setting of heated humidifier .

LCD screen

In stand by mode shows the date and hour.

Operating mode

- C for CPAP
- A for APAP

Left control button

- Allows the selection of the symbol at the lower left of the screen .
- In stand by mode gives you the access to device hour meter & summary information about the last 5 sessions.

Right control button

- Allows the selection of the symbol at the lower right of the screen .
- In stand by mode gives you access to device settings.

Ramp and selection button

- Provides access to ramp functions.
- It is also used to select the next parameter in displayed menu during the device settings.

ON/OFF button

- Turns the device ON or OFF.
- It is used to exit the device settings menu.



Setting up your DreamStar

Accessing device settings

By pressing the  button, while the device is in standby mode, to select the symbol  displayed on the screen. In the settings page display:

- The symbol  at the bottom of the screen indicates that it is impossible to change the parameters displayed.
- The symbol  at the bottom of the screen indicates that the parameters displayed can be changed:
 - Select the parameter to change by pressing the ramp button as many times as it is necessary .
 - Increase  or  decrease the parameter value by pressing the control buttons.



To return to the page header, press the ramp button  at the bottom of the page or the on/off button .

To move from one settings page to another, select the symbol  displayed on the screen by pressing the right control button . To return to the previous page, select the symbol  by pressing the left control button .

To exit the settings menu, press the on/off button  or wait 30 seconds.

Note: The following is a demonstration how to change the available settings

Humidification (if water chamber is included)

① Press  in standby mode	② The page displays: 	③ Press to select the setting 	④ The page changes: 	⑤ Increase  or decrease  the parameter value with the right or left parameter settings
--------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Settings

This page displays only if the ramp feature and (or) Pressure Alter feature and (or) the Comfort Calibration + feature were activated by the home care provider.

① DreamStar Info or Auto in CPAP mode: 	① DreamStar Auto in Auto CPAP mode: 	② Press  as many times to access the desired parameter setting	③ If CC+ is OFF, you can turn it ON by pressing the left or right parameter The press  to access the next menu parameter Increase  or decrease  the parameter value with the right or left parameter settings button
----------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Definitions

Comfort Calibration + (CC+)

If enabled by the physician or home care provider, the Comfort calibration + feature aims to increase the treatment pressure during inspiration and to decrease it during expiration to make the breathing more comfortable during the treatment.

Ramp (DreamStar Info or DreamStar Auto in CPAP mode)

If enabled by your home care provider, the ramp feature allows for a gradual rise of pressure to help you go to sleep.

Latency (DreamStar Auto in APAP mode)

If enabled by your home care provider, the latency feature allows for a gradual rise of pressure to help you go to sleep. During this time the Auto CPAP algorithm is suspended.

Pressure Alter

If enabled by your home care provider, this feature provides the option to change the prescribed pressure by ± 1 cm H₂O

Note: Pressure alter is only available in CPAP mode.

Data card

Every device is delivered with a data card. The data card allows you to download a portion of the data stored in the device's internal memory or even upload settings if your physician decides the change of the settings.

Download data

- Use the data card located at the right side of the bag.
- Insert the card carefully at the data card slot.
- Watch the progress bar at the device's screen.
- When finished remove the card.



Upload settings

- Use the card that your physician or homecare provide had programmed the settings.
- Insert the card carefully at the data card slot.
- Instantaneously the settings will be uploaded.
- When finished remove the card.



Note: The settings will be uploaded the first time you insert the data card. If you insert again the data card the device will perform compliance data download



Cleaning your DreamStar

Always refer to the device user manual for detailed instructions for the care of the device and accessories.

Daily

- Empty water reservoir and wash it with warm water and mild detergent. After the water chamber has been cleaned, you can let it soak for 15 minutes in a solution of 1 part white vinegar to 9 parts water.
- Rinse well with water to eliminate any trace of vinegar.
- Allow to dry naturally.

Note:

- *The different parts of the water reservoir may also be washed in a dishwasher (at 60°C maximum for 2 hours).*
- *To prevent micro-organisms to grow, do not allow water in the chamber to stagnate.*



Weekly

- The water reservoir as described above.
- The high efficiency filter, which is offered as an option, cannot be washed. It must be changed once a month or more frequently if visibly dirty.
- Change the filters systematically, as soon as they are torn or soiled or when the device displays a reminder to do so.



Monthly

Device

- Clean the outside of the unit using a damp tissue (cloth, paper towel) with a little water and one drop of mild detergent.
- Eliminate any trace of detergent in repeating this operation with a new tissue (cloth, paper towel) moistened with water only.
- Wipe the device thoroughly using a dry tissue (cloth, paper towel).

Air intake filter

- As described above in this page.

Water chamber

- As described above in this page.



For more information contact us at the following addresses :
customerservice@sefam-medical.com
marketing@sefam-medical.com

DreamStar™ is a trademark of SEFAM **CE** 0459

Manufacturing plant 10 Allée Pelletier Doisy - 54600 VILLERS-LES-NANCY - FRANCE
Legal entity, 144 avenue Charles de Gaulle - 92200 NEUILLY sur SEINE - FRANCE

BRO064-00/ 2012